

LENG8 Antibody



PACO10229

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LENG8

Uniprot

Q96PV6

Synonyms:

leukocyte receptor cluster (LRC) member 8; LENG8; KIAA1932; MGC40108; pp13842 ;

Immunogen:

Human LENG8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A