

LEPREL2 Antibody



PACO10234

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

P3H3

Uniprot

Q8IVL6

Synonyms:

leprecan-like 2; LEPREL2; GRCB; HSU47926; P3H3 ;

Immunogen:

Human LEPREL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A