

LIAS Antibody



PACO10247

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LIAS

Uniprot

O43766

Synonyms:

lipoic acid, synthetase; LIAS; HUSSY-01; LAS; LIP1; MGC23245 ;

Immunogen:

Human LIAS.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A