

# LILRA5 Antibody



PACO10252

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LILRA5

**Uniprot**

A6NI73

**Synonyms:**

leukocyte immunoglobulin-like receptor, subfamily A (with TM domain), member 5; LILRA5; CD85; CD85F; ILT11; LILRB7; LIR9 ;

**Immunogen:**

Human LILRA5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A