

LOH12CR1 Antibody



PACO10276

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

BORCS5

Uniprot

Q969J3

Synonyms:

loss of heterozygosity, 12, chromosomal region 1; LOH12CR1; LOH1CR12 ;

Immunogen:

Human LOH12CR1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A