

LPHN3 Antibody



PACO10282

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ADGRL3

Uniprot

Q9HAR2

Synonyms:

latrophilin 3; LPHN3; C1RL3; LEC3 ;

Immunogen:

Human LPHN3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A