

LRCH3 Antibody



PACO10285

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

LRCH3

Uniprot

Q96I18

Synonyms:

leucine-rich repeats and calponin homology (CH) domain containing 3; LRCH3; FLJ20994; FLJ43245; FLJ44976; MGC4126 ;

Immunogen:

Human LRCH3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A