

LRPPRC Antibody



PACO10289

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

LRPPRC

Uniprot

P42704

Synonyms:

leucine-rich PPR-motif containing; LRPPRC; CLONE-23970; GP130; LRP130; LSFC ;

Immunogen:

Human LRPPRC.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A