

LRRC1 Antibody



PACO10290

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LRRC1

Uniprot

Q9BTT6

Synonyms:

leucine rich repeat containing 1; LRRC1; FLJ10775; FLJ11834; LANO; dJ523E19.1 ;

Immunogen:

Human LRRC1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A