

LRRC41 Antibody



PACO10291

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LRRC41

Uniprot

Q15345

Synonyms:

leucine rich repeat containing 41; LRRC41; MGC126571; MGC126573; MUF1; PP7759; RP4-636H5.2 ;

Immunogen:

Human LRRC41.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A