

# LRRC8C Antibody



PACO10293

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LRRC8C

**Uniprot**

Q8TDW0

**Synonyms:**

leucine rich repeat containing 8 family, member C; LRRC8C; AD158; DKFZp586J1119; FAD158; MGC138551 ;

**Immunogen:**

Human LRRC8C.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A