

LSM7 Antibody



PACO10304

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LSM7

Uniprot

Q9UK45

Synonyms:

LSM7 homolog, U6 small nuclear RNA associated (S. cerevisiae); LSM7; YNL147W ;

Immunogen:

Human LSM7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A