

LSP1 Antibody



PACO10305

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LSP1

Uniprot

P33241

Synonyms:

lymphocyte-specific protein 1; LSP1; WP34; pp52 ;

Immunogen:

Human LSP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A