

LYPLAL1 Antibody



PACO10316

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LYPLAL1

Uniprot

Q5VWZ2

Synonyms:

lysophospholipase-like 1; LYPLAL1; FLJ99730; KIAA1238 ;

Immunogen:

Human LYPLAL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A