

LYZL4 Antibody



PACO10320

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LYZL4

Uniprot

Q96KX0

Synonyms:

lysozyme-like 4; LYZL4;1810009N24Rik; LYC4; MGC26768 ;

Immunogen:

Human LYZL4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A