

MAD1L1 Antibody



PACO10326

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MAD1L1

Uniprot

Q9Y6D9

Synonyms:

MAD1 mitotic arrest deficient-like 1 (yeast); MAD1L1; HsMAD1; MAD1; PIG9; TP53I9; TXBP181 ;

Immunogen:

Human MAD1L1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A