

# MAD2L2 Antibody



PACO10329

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MAD2L2

**Uniprot**

Q9UI95

**Synonyms:**

MAD2 mitotic arrest deficient-like 2 (yeast); MAD2L2; MAD2B; REV7 ;

**Immunogen:**

Human MAD2L2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A