

# MAGEA10 Antibody



PACO10334

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MAGEA10

**Uniprot**

P43363

**Synonyms:**

melanoma antigen family A, 10; MAGEA10; CT1.10; MAGE10; MGC10599 ;

**Immunogen:**

Human MAGEA10.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A