

# MAGEB18 Antibody



PACO10339

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MAGEB18

**Uniprot**

Q96M61

**Synonyms:**

melanoma antigen family B, 18; MAGEB18; MGC33889 ;

**Immunogen:**

Human MAGEB18.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A