

MAGEB2 Antibody



PACO10340

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MAGEB2

Uniprot

O15479

Synonyms:

melanoma antigen family B, 2; MAGEB2; CT3.2; DAM6; MAGE-XP-2; MGC26438 ;

Immunogen:

Human MAGEB2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A