

MAGEB4 Antibody



PACO10341

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MAGEB4

Uniprot

O15481

Synonyms:

melanoma antigen family B, 4; MAGEB4; CT3.6; MGC33144 ;

Immunogen:

Human MAGEB4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A