

MAGEC3 Antibody



PACO10343

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MAGEC3

Uniprot

Q8TD91

Synonyms:

melanoma antigen family C, 3; MAGEC3; CT7.2; HCA2; MAGE-C3; MAGEC4; MGC119270; MGC119271 ;

Immunogen:

Human MAGEC3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A