

MAGED2 Antibody



PACO10345

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

MAGED2

Uniprot

Q9UNF1

Synonyms:

melanoma antigen family D, 2; MAGED2;11B6; BCG1; HCA10; JCL-1; MAGE-D2; MAGED; MGC8386 ;

Immunogen:

Human MAGED2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A