

MAN2B2 Antibody



PACO10354

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MAN2B2

Uniprot

Q9Y2E5

Synonyms:

mannosidase, alpha, class 2B, member 2; MAN2B2; KIAA0935 ;

Immunogen:

Human MAN2B2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A