

MANBAL Antibody



PACO10355

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MANBAL

Uniprot

Q9NQG1

Synonyms:

mannosidase, beta A, lysosomal-like;

Immunogen:

Human MANBAL.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A