

MAP2K3 Antibody



PACO10362

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MAP2K3

Uniprot

P46734

Synonyms:

mitogen-activated protein kinase kinase 3; MAP2K3; MAPKK3; MEK3; MKK3; PRKMK3 ;

Immunogen:

Human MAP2K3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A