

# MAP3K10 Antibody



PACO10365

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MAP3K10

**Uniprot**

Q02779

**Synonyms:**

mitogen-activated protein kinase kinase kinase 10; MAP3K10; MLK2; MST ;

**Immunogen:**

Human MAP3K10.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A