

MAP3K11 Antibody



PACO10366

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

MAP3K11

Uniprot

Q16584

Synonyms:

mitogen-activated protein kinase kinase kinase 11; MAP3K11; MGC17114; MLK-3; MLK3; PTK1; SPRK ;

Immunogen:

Human MAP3K11.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A