

MAP4K2 Antibody



PACO10372

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MAP4K2

Uniprot

Q12851

Synonyms:

mitogen-activated protein kinase kinase kinase kinase 2; MAP4K2; BL44; GCK; RAB8IP ;

Immunogen:

Human MAP4K2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A