

MAPK8IP2 Antibody



PACO10382

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MAPK8IP2

Uniprot

Q13387

Synonyms:

mitogen-activated protein kinase 8 interacting protein 2; MAPK8IP2; IB2; JIP2; PRKM8IPL ;

Immunogen:

Human MAPK8IP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A