

MARCKSL1 Antibody



PACO10392

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MARCKSL1

Uniprot

P49006

Synonyms:

MARCKS-like 1; Marcksl1; AL022768; AW215397; AW536807; D4Bc1; F52; MacMARCKS; Macs-2; Macs-3; Macs2; Macs3; Mlp; Mrp ;

Immunogen:

Human Marcksl1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A