

MARK1 Antibody



PACO10393

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MARK1

Uniprot

Q9P0L2

Synonyms:

MAP/microtubule affinity-regulating kinase 1; MARK1; KIAA1477; MARK; MGC126512; MGC126513 ;

Immunogen:

Human MARK1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A