

MATR3 Antibody



PACO10402

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MATR3

Uniprot

P43243

Synonyms:

matrin 3; MATR3; DKFZp686K0542; DKFZp686K23100; KIAA0723; MGC9105 ;

Immunogen:

Human MATR3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A