

MCC Antibody



PACO10412

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MCC

Uniprot

P23508

Synonyms:

mutated in colorectal cancers; MCC; DKFZp762O1615; FLJ38893; FLJ46755; MCC1 ;

Immunogen:

Human MCC.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A