

MCCC2 Antibody



PACO10414

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MCCC2

Uniprot

Q9HCC0

Synonyms:

methylocrotonoyl-Coenzyme A carboxylase 2 (beta); MCCC2; MCCB ;

Immunogen:

Human MCCC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A