

MCOLN3 Antibody



PACO10426

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MCOLN3

Uniprot

Q8TDD5

Synonyms:

mucolipin 3; MCOLN3; FLJ11006; FLJ36629; MGC71509; TRP-ML3; TRPML3 ;

Immunogen:

Human MCOLN3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A