

MCTP2 Antibody



PACO10428

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MCTP2

Uniprot

Q6DN12

Synonyms:

multiple C2 domains, transmembrane 2; MCTP2; FLJ11175; FLJ33303; MGC111402 ;

Immunogen:

Human MCTP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A