

MECR Antibody



PACO10436

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MECR

Uniprot

Q9BV79

Synonyms:

mitochondrial trans-2-enoyl-CoA reductase; MECR; CGI-63; FASN2B; NRBF1 ;

Immunogen:

Human MECR.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A