

MED19 Antibody



PACO10441

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MED19

Uniprot

A0JLT2

Synonyms:

mediator complex subunit 19; MED19; DT2P1G7; LCMR1 ;

Immunogen:

Human MED19.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A