

MED26 Antibody



PACO10443

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MED26

Uniprot

O95402

Synonyms:

mediator complex subunit 26; MED26; CRSP7; CRSP70 ; ARC70;

Immunogen:

Human MED26.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A