

MED28 Antibody



PACO10444

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MED28

Uniprot

Q9H204

Synonyms:

mediator complex subunit 28; MED28;1500003D12Rik; DKFZp434N185; EG1; magicin ;

Immunogen:

Human MED28.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A