

METTL7B Antibody



PACO10468

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

METTL7B

Uniprot

Q6UX53

Synonyms:

methyltransferase like 7B; METTL7B; MGC17301 ;

Immunogen:

Human METTL7B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A