

# MFN1 Antibody



PACO10473

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MFN1

**Uniprot**

Q8IWA4

**Synonyms:**

mitofusin 1; MFN1; DKFZp762F247; FLJ20693; MGC41806; hfzo1; hfzo2 ;

**Immunogen:**

Human MFN1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A