

MICAL1 Antibody



PACO10488

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IP

Recommended dilutions:

Protein Background:

Gene ID:

MICAL1

Uniprot

Q8TDZ2

Synonyms:

microtubule associated monooxygenase, calponin and LIM domain containing 1; MICAL1; DKFZp434B1517; FLJ11937; FLJ21739; MICAL; MICAL-1; NICAL ;

Immunogen:

Human MICAL1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A