

# MIPOL1 Antibody



PACO10499

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MIPOL1

**Uniprot**

Q8TD10

**Synonyms:**

mirror-image polydactyly 1; MIPOL1; DKFZp313M2036; MGC34010 ;

**Immunogen:**

Human MIPOL1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A