

CENPU Antibody



PACO10508

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CENPU

Uniprot

Q71F23

Synonyms:

MLF1 interacting protein; MLF1IP; CENP-50; CENP-U; CENP-U(50); CENP50; CENPU; FLJ23468; KLIP1; PBIP1 ;

Immunogen:

Human MLF1IP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A