

# MMAB Antibody



PACO10519

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MMAB

**Uniprot**

Q96EY8

**Synonyms:**

methylnalonic acid, ria (cobalamin deficiency) cblB type; MMAB; ATR; MGC20496; cblB ;

**Immunogen:**

Human MMAB.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A