

MMACHC Antibody



PACO10520

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MMACHC

Uniprot

Q9Y4U1

Synonyms:

methylmalonic acid, ria (cobalamin deficiency) cblC type, with homocystinuria; MMACHC; DKFZp564I122; FLJ25671; RP11-291L19.3; cblC ;

Immunogen:

Human MMACHC.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A