

# MMP1 Antibody



PACO10521

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MMP1

**Uniprot**

P03956

**Synonyms:**

matrix metalloproteinase 1 (interstitial collagenase); MMP1; CLG; CLGN ;

**Immunogen:**

Human MMP1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A