

MOBP Antibody



PACO10539

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MOBP

Uniprot

Q13875

Synonyms:

myelin-associated oligodendrocyte basic protein; MOBP; MGC87379 ;

Immunogen:

Human MOBP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A